

ISPYB / MX CUBE JANUARY 2018

PRESNTATION OF PROCESSING AND
REPROCESSING RESULTS AT DLS

OVERVIEW

- Interactions of user with SynchWeb
- Auto-processing results
- Reprocessing - creating job, submission, results

AUTO PROCESSING - REPORT

17-01-2018 12:59:52 - 2018_01_17/rot_calib/thau/thau_1/thau_1_3_####.cbf

Sample: thau_1	Flux: 1.49e+11
Ω Start: 0.0°	Ω Osc: 0.10°
Ω Overlap: 0°	No. Images: 1200
ϕ: 45° χ: 45°	Resolution: 1.80Å
Wavelength: 0.9762Å	Exposure: 0.050s
Transmission: 9.99%	Beamsize: 80x20μm

Type: SAD

Comment: (-354,-1827,-391) (-67,-162,-272) Aperture: Large dials.align_crystalAperture: Large

Auto Processing Fast DP: ✓ Xia2: ✓✓✓ MultiXia2: ✓✓ autoPROC: ✓

Downstream Processing Fast EP: ✓ Big EP: ?? Dimple: ? MrBUMP: ✘

0 20 40 60 80 100
1200 1100 1000 900 800 700 600
1.35 1.40 1.45 1.50 1.55 1.60 1.65 1.70

Spots Bragg Res

17-01-2018 12:57:58 - 2018_01_17/rot_calib/thau/thau_1/thau_1_2_####.cbf

Sample: thau_1	Flux: 1.50e+11
Ω Start: 0.0°	Ω Osc: 0.10°
Ω Overlap: 0°	No. Images: 1200
ϕ: 0° χ: 45°	Resolution: 1.80Å
Wavelength: 0.9762Å	Exposure: 0.050s
Transmission: 9.99%	Beamsize: 80x20μm

Type: SAD


Comment: (-606,-1604,-483) (-67,-162,-272) Aperture: Large dials.align_crystalAperture: Large

Auto Processing Fast DP: ✓ Xia2: ✓✓✓ MultiXia2: ✓✓ autoPROC: ✓


Downstream Processing Fast EP: ✓ Big EP: ?? Dimple: ? MrBUMP: ✘

0 20 40 60 80 100
1300 1200 1100 1000 900 800 700 600
1.35 1.40 1.45 1.50 1.55 1.60 1.65 1.70


Spots Bragg Res



DETAILS




DETAILS




REPROCESSING



REPROCESSING






CREATE REPROCESSING JOB

- 1 entry in ProcessingJob (describing the job)
- 1 or more entries in ProcessingJobImageSweep
(selecting data to be used for the job)
- 0 or more entries in ProcessingJobParameter (arbitrary key/value pairs)

in the database, plus sending a trigger message to ActiveMQ which references the ProcessingJob row.



TABLES

ProcessingJobImageSweep

processingJobImageSweepId - Autoincrementing unique ID
processingJobId - Job this entry belongs to
dataCollectionId - because although a ProcessingJob is already connected to one DCID, data may come from other DCIDs as well

startImage

endImage

ProcessingJobParameter

processingJobParameterId - Autoincrementing unique ID
processingJobId - Job this entry belongs to
parameterKey
parameterValue



```
{  
  "environment": {  
    "ID": "38c1219c-8ae9-4b22-8e62-8fd5e1594194",  
    "ispyb_autoprocprogram_id": 68626954  
  },  
  "recipe": {  
    "1": {  
      "output": [  
        2  
      ],  
      "parameters": {  
        "cmdline": "xia2 via zocalo",  
        "environment": {},  
        "ispyb_command": "register_processing",  
        "program": "xia2-dials",  
        "rpid": "4024",  
        "store_result": "ispyb_autoprocprogram_id"  
      },  
      "queue": "ispyb_connector",  
      "service": "DLS ISPyB connector"  
    },  
    "2": {  
      "job_parameters": {  
        "dcid": "2227029",  
        "ispyb_parameters": {},  
        "results_directory": "/dls/i03/data/2018/cm19644-1/processed/2018_01_17/  
rot_calib/thau/thau_1/thau_1_2_/86449557-7045-44dd-8892-0ccff51815e4/ispyb-xia2-  
dials",  
        "timeout": null,  
        "working_directory": "/dls/i03/data/2018/cm19644-1/tmp/zocalo/2018_01_17/  
rot_calib/thau/thau_1/thau_1_2_/86449557-7045-44dd-8892-0ccff51815e4/ispyb-xia2-  
dials",  
      }  
    }  
  }  
}
```



```
xia2 pipeline=dials min_images=3 failover=True \
project=AUTOMATIC crystal=DEFAULT read_all_image_headers=False atom=S \
image=/dls/.../rot_calib/thau/thau_1/thau_1_2_0801.cbf:801:1200 \
image=/dls/.../rot_calib/thau/thau_1/thau_1_2_0001.cbf:1:400 \
image=/dls/.../rot_calib/thau/thau_1/thau_1_2_0401.cbf:401:800
```

RESULTS

17-01-2018 12:57:58 - ▶ 2018_01_17/rot_calib/thau/thau_1/thau_1_2_####.cbf

Sample: thau_1	Flux: 1.50e+11
Ω Start: 0.0°	Ω Osc: 0.10°
Ω Overlap: 0°	No. Images: 1200
ϕ: 0° χ: 45°	Resolution: 1.80Å
Wavelength: 0.9762Å	Exposure: 0.050s
Transmission: 9.99%	Beamsize: 80x20μm
Type: SAD	

Comment: (-606,-1604,-483) (-67,-162,-272) Aperture:
Large dials.align_crystalAperture: Large

Auto Processing

Fast DP: ✓ Xia2: ✓ ✓ MultiXia2: ✓ ✓ autoPROC: ✓

Fast DP XIA2 3d XIA2 DIALS XIA2 3dii autoPROC XIA2 3dii XIA2 DIALS xia2 via zocalo

Beam Centre	X	Y
Start	211.42	206.11
Refined	211.46	206.04
Δ	-0.04	0.07

Space Group	A	B	C	α	β	γ
P 41 21 2	57.82	57.82	150.55	90.00	90.00	90.00

Shell	Observations	Unique	Resolution	Rmeas	I/sig(I)	CC Half	Completeness	Multiplicity	Anom Completeness	Anom Multiplicity	CC Anom
outerShell	703	512	1.36 - 1.38	0.348	1.4	0.8	18.7	1.4	2.7	1.2	0.0
innerShell	46837	3116	3.69 - 57.87	0.048	38.4	1.0	100.0	15.0	100.0	8.9	-0.2
overall	572636	49004	1.36 - 57.82	0.076	17.2	1.0	87.8	11.7	80.9	6.5	-0.3

Downstream Processing

Fast EP: ✓ Big EP: ? ? Dimple: ? MrBUMP: ✗